




ANTIBACTERIAL DRESSING

MEDIHONEY®

Antibacterial Honey Tulle 3-Ply Dressing



FEATURES

- Antibacterial *Leptospermum* (Manuka) Honey (minimum 20g)
- Non adherent

INDICATIONS

- Leg/foot ulcers
- Pressure ulcers
- Infected wounds
- Sloughy wounds
- Necrotic wounds
- Malodorous wounds
- Donor and recipient graft sites
- Burns
- Surgical wounds
- Diabetic wounds
- Abrasions

SUPERIOR WOUND BED PREPARATION

SUPERIOR WOUND BED PREPARATION

MEDIHONEY® Antibacterial Honey Tulle Dressing protects the wound by creating a barrier against wound pathogens, including antibiotic resistant strains, and therefore reducing the risk of infection. The osmotic action produces an outflow of body fluid which assists the removal of wound bacteria, endotoxins, debris and slough, providing a cleaner wound, rapidly removing malodour and helping to reduce the inflammatory response, oedema and exudate levels. Granulation and epithelialisation are enhanced through provision of the optimal healing environment.

MEDIHONEY® Antibacterial Honey Tulle Dressing


creates a microenvironment that supports healing and provides an antibacterial barrier.



ANTIBACTERIAL DRESSING

MEDIHONEY®

Antibacterial Honey Tulle 3-Ply Dressing



DIRECTIONS FOR USE

1. If required, protect the wound edges with MEDIHONEY® Barrier Cream.
2. Place the dressing onto the wound surface ensuring that the dressing is in contact with the wound bed. Dressings can be placed side by side to cover large areas or for deep wounds multiple layers may be used. Alternatively the wound may be filled with MEDIHONEY® Antibacterial Medical Honey® or MEDIHONEY® Antibacterial Wound Gel® before applying MEDIHONEY® Antibacterial Honey Tulle Dressing.
3. A suitable sterile pad or dressing should be applied to cover the MEDIHONEY® Antibacterial Honey Tulle Dressing. The dressing should be sufficiently absorbent to manage wound exudate.
4. The wound may be cleaned between dressing changes in accordance with standard wound cleaning protocol.

Dressings may be folded or cut with sterile scissors to the shape of the wound.

Frequency of change: MEDIHONEY® Antibacterial Honey Tulle Dressing requires changing when the dressing has been diluted by wound exudate. The dressing may be left in place for up to 7 days depending on wound exudate. Heavily exudating wounds require close monitoring and may require daily dressing changes.

ITEM#	PRODUCT NAME	PIP CODE	NHS CODE
796	10cm x 10cm, Tulle Dressing	346-3346	EJE067