


Advancing the **GOLD STANDARD** of care to the next generation.



**QUICKER APPLICATION,
GREATER COMFORT.**

TCC-EZ® lightweight, customized boots simplify the process and provide additional support and stability.

TCC-EZ® Boot Dimensions

Measurements	TCC-EZ® Regular Boot TCC2-1100	TCC-EZ® Large Boot* TCC2-1114	TCC-EZ® Extra Large Boot* TCC2-1124
Length	28.9 cm (11.4 in)	31 cm (12.25 in)	34.9 cm (13.75 in)
Width at Ankle	10 cm (3.94 in)	10.5 cm (4.11 in)	10.9 cm (4.29 in)
Width at Toe	11.8 cm (4.66 in)	12 cm (4.76 in)	13 cm (5.14 in)
Approximate Shoe Size			
Men's	6 - 10	10.5 - 13 wide	13.5 - 18 X-wide
Women's	5 - 11.5	11.5+ wide	11.5+ X-wide
Patient Weight	Up to 125 kg (275 lbs)	125 - 181 kg (275-400 lbs)	Up to 181 kg (400 lbs)

*Patients who wear wide shoes may need to be evaluated for a Large or Extra Large Boot regardless of foot length.



TCC-EZ® offers Large and Extra Large Charcot Boots with innovative cushion insert specially designed to accommodate patients who have developed rocker bottom deformities.

TCC-EZ® Charcot Boot Dimensions

Measurements	TCC-EZ® Large Charcot TCC2-1116	TCC-EZ® Extra Large Charcot TCC2-1126
Length	31 cm (12.25 in)	34.9 cm (13.75 in)
Width at Ankle	10.5 cm (4.11 in)	10.9 cm (4.29 in)
Width at Toe	12 cm (4.76 in)	13 cm (5.14 in)
Approximate Shoe Size		
Men's	14 wide	14 - 18 X-wide
Women's	14 wide	14+ X-wide
Patient Weight	Up to 181 kg (400 lbs)	Up to 181 kg (400 lbs)



104 Shorting Road, Toronto, ON M1S 3S4 Canada
P: 416.299.4003 | Toll Free: 800.387.5302
E: customerservice@dermasciences.com

www.TCCEZ.com
www.TCCPatient.com

A TISSUE REGENERATION COMPANY

Refer to the complete instructions for use as well as boot instructions for use. Never walk without the boot.
1. Armstrong DG, et al. Off-loading the diabetic foot wound. Diabetes Care 24:1019-1022, 2001.
2. Bloomgarden ZT: American Diabetes Association 60th Scientific Sessions, 2000. Diabetes Care 24:946-951, 2001.
© 2014 Derma Sciences, Inc. All rights reserved. TCC Boot-001