

ANTIMICROBIAL MEDIHONEY[®] IN ABSORBENT HYDROGEL



ABSORBENT

63%
OF THIS PRODUCT
IS 100% ACTIVE
MANUKA HONEY

*Active Leptospermum Honey

OVERVIEW

Medihoney[®] HCS is an all in one dressing combining the unique properties of medical grade Manuka honey with a hydrogel sheet containing superabsorbent polymers.

FEATURES AND BENEFITS

- 63% active Manuka[®] honey offers proven antimicrobial and anti-inflammatory actions
- An all in one low profile dressing (adhesive version) ensures ease of application, suitable under compression
- Gentle release of honey into the wound provides a low pH environment and autolytic debridement, encouraging faster wound healing
- Absorption of exudate into the dressing is as a solid gel, thus maximises wear time and minimising dressing changes
- Highly flexible, transparent dressing offers a soothing and cooling effect on application, improving patient comfort
- Atraumatic removal minimises disruption to the wound, reducing pain
- Can be cut to treat different shaped wounds and around exit tube sites

INDICATIONS

- Diabetic Foot Ulcers
- Leg Ulcers (venous, arterial and mixed aetiology ulcers)
- Pressure Ulcers
- 1st and 2nd degree burns
- Donor sites
- Traumatic Wounds
- Surgical and post op wounds
- Skin tears
- Suitable for neonatal and paediatric patients



ANTIMICROBIAL MEDIHONEY[®] IN ABSORBENT HYDROGEL



| Code | Description | Product Dimensions | Qty/Box | PIP Code | NHS Code |
|------|--|--|---------|----------|----------|
| 780 | Medihoney [®] HCS Non-Adhesive | 6cm x 6cm | 10 | 382-2491 | ELM200 |
| 781 | Medihoney [®] HCS Non-Adhesive | 11cm x 11cm | 10 | 382-2525 | ELM201 |
| 784 | Medihoney [®] HCS Non-Adhesive | 20cm x 20cm | 5 | | EKB090 |
| 785 | Medihoney [®] HCS Non-Adhesive | 20cm x 30cm | 2 | | EKB089 |
| 782 | Medihoney [®] HCS Adhesive | Inner 7.2cm x 7.2cm, outer 11cm x 11cm | 10 | 382-2509 | ELM202 |
| 783 | Medihoney [®] HCS Adhesive | Inner 11.5cm x 11.5cm, outer 15cm x 15cm | 10 | 382-2517 | ELM203 |
| 787 | Medihoney [®] HCS Adhesive Surgical | Inner 4.5cm x 16.5cm, outer 7.5cm x 20cm | 10 | | EKB094 |

Paediatric Facial Lesions

A 15-year-old female patient developed severe right-sided facial cellulitis with vesicular and ulcerative lesions measuring 8cm x 8cm with a depth ranging from 0.1cm to 0.3cm. The wound bed comprised 100% slough or hyperpigmented eschar. There were multiple full-thickness and partial-thickness open wounds scattered over the right side of face with oedema and erythema. Tests revealed no clear aetiology to explain the infection.



Figure 1:
Before application



Figure 2:
At start of therapy with
Medihoney[®] HCS



Figure 3:
4 days after the start of therapy.
All eschar removed with 50% slough
remaining and no residual crusting



Figure 4:
4 weeks after start of therapy, all
slough removed and the wounds
were almost completely granulated
or epithelialised